

## PSU-2T, Mini-shaker for immunology

Shaker **PSU-2T** provides adjustable mixing of reagents in microtest plates. The shaker is a compact and user-friendly device; it occupies small space on a desk. The instrument is designed for immunoassay. The shaker is ideal for personal use.

Shaker can be used in cold rooms or incubators, operating at ambient temperature range +5°C to +40°C. Low voltage external power supply (12 V) provides electrical safety in humid environment.

### Specifications:

Speed range	150–1000 rpm
Digital timer	1 min–23 hrs 59 min
Orbit diameter	2 mm
Overall dimensions	220 x 205 x 90 mm
Weight, not more	2.5 kg
External power supply	DC 12 V, 500 mA


### Platforms for microtest plates:

**A** IPP-2 (standard platform) for 2 microtest plates

**B** IPP-4 (optional platform) for 4 microtest plates

### Catalogue number:

<b>PSU-2T</b> with platform <b>IPP-2</b>	BS-010102-AAG
<b>IPP-4</b>	BS-010102-AK



Product video is available on the website

### **A** Platform IPP-4



### **B** Platform IPP-2

